

Computer Assisted Neural Probe Insertion

Shizra Tariq
04/10/2025



UNIVERSITY OF MINNESOTA
Driven to DiscoverSM

Introduction


Goal

- Automate probe insertion into mouse brain
- Use computer vision for real-time guidance
- Increase accuracy and safety
- Reduce human error

Why we need to solve this problem???


Challenges

- Manual insertion is slow and skill-dependent
- Small mistakes → probe bending/breaking or tissue damage
- Difficult to visualize and control probe position
- Need for consistent and precise insertions




Available Datasets

Two kinds of datasets that were available when I joined the project:



Dataset 1: zoomed-in visual setup

Probes are more visible but hard to detect the insertion point due to lack of surrounding reference



Dataset 2: Wider Camera View

probes are comparatively less visible but provides partial view of the headplate

Probe Detection


- Our first goal was to detect the probe during insertion.
- We focused on identifying both the probe shaft and probe tip using camera views.
- Once the probe was detected, the next step was to classify whether it was straight or bent.
- Detecting bending is important because it means the probe is not aligned correctly or hitting resistance.
- Early detection of bending helps prevent probe damage or breakage.

Proposed Solution


Designed a two-stage pipeline:

- Segmentation
- Classification



Methodology




Results



Detection accuracy for bent probes was comparatively lower due to dataset imbalance






New Dataset




- Collected new dataset after camera change
- Increased noise and background complexity
- Change enabled clear visibility of headplate
- Headplate provides reference for insertion point estimation

Results



Insertion Point Estimation

- Goal: locate insertion points for automated probe placement.
- Manual insertion requires precise identification of brain coordinates.
- Automation computes points relative to stable anatomical landmarks.
- Requires a stable reference for accurate calculation.
- Accurate detection is critical to avoid tissue damage.




Initial Approach

Designed reference-based localization:

- **Headplate Segmentation Approach:**

Segment headplate → use as fixed reference → calculate insertion coordinates



- **Alternative Approach:**

Detect [anatomical features/markers] → establish reference frame → compute insertion points

Why it failed?


Initial Assumption: Camera would remain fixed throughout all trials

- This would allow a one-time calibration using the headplate reference



Reality: Camera position is only static **per trial**, not across trials

- Camera repositioning between trials changes the reference frame
- Each new trial requires recalibration relative to the headplate


- Camera repositioning alters headplate position, scale, and orientation, invalidating calibrations.
- Small segmentation errors (2–3 px) cause large deviations; sub-millimeter accuracy is required.
- Converting 2D to 3D coordinates on the curved headplate requires precise per-camera calibration.
- Lack of fiducial markers makes repeatable reference point identification challenging.




Exo Setup






Why 3D conversion is necessary?






Camera Calibration



Zhang's calibration method



Detection with ArUco tags



Find marker
corners and tag
ID

Convert to 3D
camera coordinates

Estimate Marker Poses

$t = \begin{bmatrix} t_x \\ t_y \\ t_z \end{bmatrix}$ $rvec = \begin{bmatrix} r_x \\ r_y \\ r_z \end{bmatrix}$

Insertion Point Calibration



Process

1. Manually identify 4 insertion points on the camera view while all ArUco markers are visible
2. converts 2D pixel coordinates to 3D headplate coordinates using the ArUco marker reference frame




Reference Coordinate System


- **Origin:** Centroid of 4 ArUco tags (44, 32, 29, 49)
- **X-axis:** Direction from tag 44 → 32 (projected to plane)
- **Z-axis:** Average normal vector from all 4 tags
- **Y-axis:** Cross product completing right-handed system

Output

#2 Detection with Aruco tags




Calibrating Insertion points




Live Detection

#2 Detection with 4 Aruco tags



Calibrating Insertion points



Live Detection

Trajectory Planning

The Goal: Go From Point A to Point B

Waypoint Structure

position: 3D coordinates (x, y, z)
orientation: Approach direction vector
speed: Velocity scaling factor
waypoint_type: Stage identifier

Trajectory Components

1. Waypoint List: Ordered sequence of positions
2. Distance Metrics: Total path length calculation
3. Time Estimation: Speed-based duration prediction

Waypoint 1
Current Probe
position

Point A = The 3D
position of the probe tip

Direction Vector: $V_seg1 = P_W2 - P_W1$
Distance (Euclidean): $D_seg1 = ||V_seg1||$

Point B = The 3D position
of the target

Waypoint 3
Insertion depth

$D_total = D_seg1 + D_seg2$

$$D = \sqrt{(X_W1 - X_W2)^2 + (Y_W1 - Y_W2)^2 + (Z_W1 - Z_W2)^2}$$

Trajectory Planning

The Goal: Go From Point A to Point B

Point A = The 3D position
of the probe tip.

Point B = The 3D position
of the target

Waypoint 1
Current Probe
position

Waypoint 2
Target surface


Waypoint 3
Insertion depth

Direction Vector: $V_{seg1} = P_{W2} - P_{W1}$
Distance (Euclidean): $D_{seg1} = ||V_{seg1}||$

$D_{total} = D_{seg1} + D_{seg2}$

$$D = \sqrt{(X_{W1} - X_{W2})^2 + (Y_{W1} - Y_{W2})^2 + (Z_{W1} - Z_{W2})^2}$$

Results



Trajectory Visualization:

- Completed waypoints (gray)
- Current waypoint (cyan)
- Approach phase (yellow)
- Insertion phase (red)
- Retraction phase (purple)

Next Steps

- Validate Accuracy: Perform trial insertions to quantify system accuracy and repeatability.
- Automate GUI: Integrate vision software with the New Scale motor controller.
- Add Multi-Camera Tracking: Implement a stereoscopic camera system to improve 3D tracking.
- Scale to Four Probes: Expand the framework to manage four automated insertions simultaneously.

Impact

Higher-Quality, Reproducible Data

- Computer vision enables consistent and precise insertions.
- Ensures probes are in the correct location, leading to reliable scientific data.

Increased Safety & Cost-Effectiveness

- Real-time feedback (like bend detection) prevents probe breakage.
- Saves money and minimizes tissue damage, improving animal welfare.

Accelerated Scientific Discovery

- Reduces human error and speeds up the entire experimental process.
- Opens the door to new, complex experiments (e.g., multi-probe insertions).

Special Thanks

Travis and **Suahasa** for helping and guiding me throughout the project

Michael for explaining and entertaining all of my questions to understand the previous work and the setup

Thank you

Questions?



UNIVERSITY OF MINNESOTA

Driven to DiscoverSM